Planopol-3

Instruction Manual
Gebrauchsanweisung
Mode d'emploi

Serial No.:

Voltage:

*1991.09

In technical enquiries and in orders for spare parts, please state apparatus No. and voltage (see cover).

All information and data in this Instruction Manual were valid when it was printed. The Struers policy being, however, to make current improvements of our products, we reserve ourselves the right to modify our products without notice. Besides, we ask you to note that this Instruction Manual concerns apparatus and all accessories.

Consequently, it may happen that the Manual mentions equipment which is not of relevance for your purpose at the present time.

- 1. Description**
- 2. Technical data**
- 3. Accessories**
- 4. Contents of packing case**
- 5. Installation**
 - 5.1 Setting up**
 - 5.2 Connection to mains**
 - 5.3 Connection to water**
 - 5.4 Mounting of Conversion Kit**
- 6. Operation**
- 7. Maintenance**
- 8. Trouble-shooting**
- 9. Spare parts (Appendix)**
- 10. Illustrations (Appendix)**

Struers A/S

Valhøjs Allé 176

DK-2610 Rødovre/Copenhagen

Denmark

☎ Telephone +45 31 70 80 90

☎ Fax +45 31 41 65 44

Telex 19625

1. Description

Planopol-3 is a standard machine for lapping, wet grinding and diamond and alumina polishing of metallographic, ceramic and mineralogical specimens.

It has a circular turntable accepting grinding, lapping and/or polishing discs of diameters up to 300 mm.

Supported by precision ball bearings, the turntable is driven via a V-belt drive by a powerful two-speed electric motor, which ensures smooth and noiseless running.

Motor and turntable are mounted on a steady steel frame also supporting the cabinet. The latter is provided with a needle valve for adjusting water admission, a flexible sprinkler tube and a spill pan with drain.

On the control panel are arranged, to the right, three push buttons with light indication for setting the speed (high/low) and for stop.

The mains switch is easily accessible to the right on the front of the cabinet.

2. Technical data

<i>Voltage/frequency</i>	<i>1-phase machines</i> 115 V, 60 Hz 220 V, 50 Hz
	<i>3-phase machines</i> 3 x 220 V, 50 Hz 3 x 380 V, 50 Hz 3 x 415 V, 50 Hz 3 x 440 V, 50 Hz 3 x 500 V, 50 Hz 3 x 220 V, 50 Hz 3 x 440 V, 60 Hz 3 x 480 V, 60 Hz
<i>Motor output</i>	<i>1-phase</i> 50 Hz 370, 220 W 60 Hz 450, 300 W
	<i>3-phase</i> 50 Hz 550, 300 W 60 Hz 650, 370 W
<i>Rotational speed</i>	<i>Standard</i> 150, 300 rpm
	<i>Conversion Kit</i> 75, 150 rpm or 300, 600 rpm or 600, 1200 rpm
<i>Polishing/grinding disc</i>	Max. 300 mm
<i>Dimensions</i>	Height: 340 mm
	Width: 462 mm
	Depth: 743 mm
<i>Weight</i>	53 kg
<i>Safety standard</i>	IEC 204/ (VDE 0113) EN 60204-1

3. Accessories

A complete range of different discs of diameters between 200 and 300 mm is available, including discs for grinding with grinding paper, lapping on cast iron with a variety of surfaces, fine lapping on plastics or special lapping discs, as well as polishing discs of PVC and aluminum. In addition, a complete program of consumables are available, such as polishing cloths, diamond paste or spray, SiC paper or powder, alumina, etc. Please refer to separate leaflets.

Planopol-3 has been prepared for mounting of PdM-Force (fig. 2), a motor-driven sample mover for simultaneous lapping and/or polishing of 1-9 specimens.

For automatic preparation of from 3 to 12 specimens simultaneously also Pedemax-2 (fig. 3) may be mounted on Planopol-3 (3-phase model), thus making it possible to grind, lap and polish a number of specimens at a time. Both PdM-Force and Pedemax-2 reduce the preparation time and provide reproducible results. For details you are referred to separate leaflets.

Conversion Kit

Planopol-3 is supplied as standard with 150/300 rpm. By mounting of a Conversion Kit (codeword: PLATT) the following rotational speeds may be obtained: 75/150, 300/600, and 600/1200 rpm.

4. Contents of packing case

- 1 Planopol-3
- 1 Protection ring
- 1 Drain elbow pipe
- 1 Drain hose
- 2 Hose clips 20-32 mm
- 1 Hose connector
- 1 Union nut for hose connector
- 1 Sealing ring for hose connector
- 1 Hose clip 13 mm
- 1 Hexagon spanner 6 mm for turntable
- 1 Fuse 0.5 AT 220 V

Conversion Kit

- 1 V-belt pulley with two grooves (fig. 5.3)
- 1 Fitted bolt (fig. 5.2)
- 1 Washer (fig. 5.4)
- 1 Adjusting plate (fig. 5.1)
- 1 V-belt Z 17 (10 x 425 mm) (60 Hz)
- 1 V-belt Z 18 (10 x 450 mm) (50 Hz)
- 1 V-belt HY-T SPZ (10 x 1312 mm) (60 Hz)
- 1 Hexagon spanner (8 mm) for fitted bolt
- 1 Hexagon spanner (2.5 mm) for V-belt pulleys
- 3 Plastic plates stating the speeds 75, 600 and 1200 rpm

5. Installation

5.1 Setting up

Planopol-3 is detached from the bottom of the packing case by removing the four screws from below, after which Planopol-3 can be placed on the bench.

5.2 Connection to mains

Before the apparatus is connected it must be checked that the voltage which is indicated on the label of the machine, is in accordance with the local voltage.

Mount a plug on the electric cable, connecting the leads as follows:

1-phase machine

Yellow/green lead to earth
Brown lead to phase
Blue lead to neutral

3-phase machine

Yellow/green lead to earth
Brown lead to phase
Black lead to phase
Black lead to phase

The turntable should rotate anti-clock wise. If this is not the case after the connection, change the direction by changing over two of the phase leads.

5.3 Connection to water

Connect the thin reinforced hose to a water outlet and fit the 13 mm hose clip. The accompanying hose connector with union nut and sealing ring may, where not already provided, be mounted on the water outlet (1/2" pipe thread, BS 2779: G 1/2).

If a recirculating cooling unit is used, connect the hose to the discharge branch of the pump and fit a hose clip.

Lead the thick hose to the drain and be very careful to place it with a steady downward slope to prevent the water from being discharged too slowly which may cause overflow or blocking of the water outlet. In order to avoid elbows on the hose, the enclosed elbow pipe can be mounted on the tubing in an appropriate place.

5.4 Mounting of Conversion Kit

The rotary speed is changed by changing the transmission ratio of motor and turntable. The V-belt pulley of the motor is provided with two grooves of different diameters. In the standard version of 150/300 rpm the drive is effected directly from the small diameter groove to the V-belt pulley of the turntable. If the motor V-belt pulley is inverted so that the drive is established from the large diameter groove, the speed of the turntable is doubled. Like the motor V-belt pulley, the extra V-belt pulley has two grooves of different diameters. By inserting this extra pulley as an intermediate wheel a speed either twice or half the one mentioned above may be obtained, depending on which way the intermediate wheel and the motor V-belt pulley are fitted.

The various combinations possible will appear from fig. 4, which also shows which V-belts to use. Always make sure that the V-belt pulleys are mounted in such a way that the belt pulls between two grooves running in the same plane.

Change of rotational speed is done as follows (after disconnection of power plug): Take out the 6 screws securing the rear section of the cabinet and lift this section off. Remove the nuts securing the motor bridge, after which the V-belt pulley may be mounted according to fig. 5. Fit the V-belt pulley as shown in fig. 4, which indicates the correct positioning of the motor V-belt pulley. The motor V-belt pulley is fastened by means of two pointed screws in the bottom of the V-belt grooves. When these screws are slackened, the V-belt pulley can be pulled off.

For exchange of the long V-belt the turntable must be taken out. This is done after removal of three screws at the top of the turntable.

Position the belts as shown in fig. 4. First tighten the short V-belt by means of the adjusting plate (fig. 5.1), then clamp the latter. Remount the motor bridge and tighten the long belt by pushing the motor rearward and clamp it into position.

It is important to ensure adequate belt tension since undertensioning will cause the belt to slip on the pulleys, impairing the drive transmitted to the turntable, whereas overtensioning produces noise.

The long belt is correctly tensioned when a force of 5 N applied halfway between the pulleys causes a bend of about 6 mm. When applied to the short belt, the same force should produce a bend of 2-3 mm.

The push button plates indicating the new speeds are exchanged by first unscrewing the metal ring round the buttons and then pressing off the green plastic plate. Fit the new indications and assemble the buttons in reverse order. Finally, remount the rear section of the cabinet.

6. Operation

Planopol-3 is started by first operating the mains switch (fig. 1.1). This actuates the red lamp in the uppermost of the three push buttons at the right side (fig. 1.2). The button in the middle of the three (fig. 1.3) starts the motor at low speed, whereas the lowest button (fig. 1.4) switches to the higher speed. Activation of these buttons causes the lamp built into each button to glow. The machine is stopped by pressing the uppermost, red push button (fig. 1.2). The push buttons may be activated in any order. When stopped by the mains switch or by a power failure, the motor does not start on reconnection of power, but has to be started by operation of the push buttons. This prevents unintentional motor start during exchange of grinding/polishing disc.

In the case of prolonged overloading of the motor a built-in thermal relay cuts out the motor, re-start being possible only after a period of cooling.

The grinding, lapping or polishing disc is mounted on the turntable by inserting the three pins of the disc in the corresponding holes in the turntable. Make sure that the contact surfaces are clean. The water flow is controlled by a needle valve (fig. 1.5), and the distribution of water is controlled by means of the flexible sprinkler tube (fig. 1.6).

Below are stated the appropriate speeds for the various preparation processes:

<i>75 rpm</i>	Lapping on cast iron discs and PAD-K pellow cloth. Diamond polishing on Petrodisc and PA-W, PAN-W pellow cloth.
<i>150 rpm</i>	Diamond polishing (grinding/lapping on Petrodisc-M).
<i>300 rpm</i>	Wet grinding on grinding paper/wheel. Grinding on Petrodisc-M or DP-Net.
<i>600 rpm</i>	Alumina polishing (automatic preparation).
<i>1200 rpm</i>	Alumina polishing (manual preparation).

The grinding paper may be affixed to a plane polishing disc by means of Spray Adhesive, or it may be placed on a special disc according to the Knuth-Rotor principle. By the latter method the disc is filled with water while at rest, after which paper and retaining ring are put on. Then the disc is made to rotate 300 rpm while the entire surface of the grinding paper is pressed gently against the disc with a specimen, or the like.

7. Maintenance

Planopol-3 is designed and constructed for many years of trouble free running. It will not need much attention, as its bearings are sealed and require no further lubrication. So maintenance will generally be limited to cleaning the exterior of the machine, the spill pan and the drain hose. Where the machine is used, for instance, for alumina polishing or SiC lapping, the spill pan should be cleaned sufficiently often to prevent drying out of the suspension and thus any risk of blocking of the hose. In order to avoid erratic disc movements, the contact surfaces between disc and turntable should be kept clean.

Now and again the V-belt may need tightening, the transmission being less efficient with inadequate belt tension. To tighten the belt, remove the rear part of the cabinet, slacken the nuts securing the motor plate and push plate and motor towards the rear until the belt tension is correct. The belt is tightened correctly when a force of 5 N, applied to the belt in the middle between the pulleys, causes a bend of 6 mm.

Exchange of the belt is effected by additionally removing the turntable and loosening the nuts securing the motor plate so as to slacken the belt, after which the belt can be taken off the turntable pulley from above. Mounting of the new belt is done in reverse order.

8. Trouble-shooting

Indication	Explanation/cause	Action
Motor does not run. Red light in stop button.	Motor overload. Thermal relay has cut out power.	Wait 5 minutes.
Motor does not run.	Internal transformer fuse blown.	Disconnect power supply. Remove rear cabinet section and replace transformer fuse.
Motor does not run, but buzzes.	One phase lacking (3-phase machine).	Check all fuses of local wiring.
Motor will only run at one speed.		Call service man.
No light in push button, otherwise normal machine function.	Bulb burnt out.	Disconnect power supply. Exchange bulb by unscrewing ring round push button. Insert new bulb and turn left.
Poor drive.	Slipping V-belt(s).	Tighten, see section 7.
Water not drained.	Hose pressed flat. Dirt trapped in hose. Presence of air stops discharge.	Adjust hose position. Clean hose. Discharge hose should have a steady downward slope.
Machine produces noise.	Grinding/polishing disc vibrates. Bearings of motor or turntable defective.	Clean contact surfaces between disc and turntable. Exchange bearings/motor.

If other problems arise, or if the above failure corrections do not help, please call for a service technician.

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Fig./ Pos.	Description	Spare part No.
<i>Electric parts</i>		
6.29	Motor 3 x 220 V/50 Hz	414 MP 29
	Motor 3 x 380 V/50 Hz	414 MP 30
	+ 3 x 440 V, 60 Hz	
	Motor 3 x 415 V, 50 Hz	414 MP 31
	Motor 3 x 440 V, 50 Hz	3940238
	Motor 3 x 500 V, 50 Hz	414 MP 32
	Motor 3 x 220 V, 60 Hz	414 MP 33
	Motor 3 x 480 V, 60 Hz	3620058
	Motor 1 x 115 V, 60 Hz	414 MP 34
	Capacitor 115 V, 60 Hz	414 MP 35
	Motor 1 x 220 V, 50 Hz	414 MP 36
	Capacitor 220 V, 50 Hz	414 MP 35
6.38	Main switch 25000100	414 MP 38
6.39	Contacteur 50 Hz	414 MP 39
	Contacteur 60 Hz	414 MP 40
6.41	Thermorelay 3 x 220 V, 50 Hz	414 MP 41
	Thermorelay 3 x 220 V, 60 Hz	414 MP 42
	Thermorelay 3 x 380 V, 50 Hz	414 MP 43
	Thermorelay 3 x 415 V, 50 Hz	414 MP 43
	Thermorelay 3 x 440 V, 60 Hz	414 MP 43
	Thermorelay 3 x 500, 550 V, 50 Hz	414 MP 46
	Thermorelay 1 x 220 V, 50 Hz	414 MP 42
	Thermorelay 1 x 115 V, 60 Hz	414 MP 48
6.49	Transformer	414 MP 49
6.50	Fuseholder	414 MP 50
6.51	Pushbutton without change-over switch, 0-red	414 MP 51
	Pushbutton without change-over switch, 150-green	414 MP 52
	Pushbutton without change-over switch, 300-green	414 MP 53
6.54	Change-over switch	414 MP 54
6.55	Bulb 24 V, 0.08 A	414 MP 55
	Fuse 0.5 AT	414 MP 56
6.57	Screwed connection for cable PG 11	414 MP 57
<i>Water inlet - outlet</i>		
6.63	Reinforced tube 5/16" x 2 m with 1/4 nipple	414 MP 63
6.58	Needle valve for water	414 MP 58
6.59	Reinforced tube 3/8" x 300 mm	414 MP 59
6.60	Adjustable metal tube 3/8" x 320 mm	414 MP 60
6.61	Outlet tube dia. 25 x 350 mm	414 MP 61
6.62	Outlet tube branch	414 MP 62

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